



RECOMMENDATIONS FOR RECOMMISSIONING FULL (MAINS) OR SINGLE SYSTEM DISINFECTION (DC / HWS)

STANDARD RECOMMISSIONING PROCEDURE (WITH LAB CULTURE ONLY)



ACCELERATED RECOMMISSIONING PROCEDURE WITH HYDROSENSE



***Always perform normal safety checks before reopening.**

Hydrosense accelerates recommissioning so you can get back to 'business as usual' quickly, and safely.

Results in 25 minutes vs 2 weeks facilitates ultra-fast validation of flushing / disinfection.

Lockdown has created many low or no maintenance systems with low flow. Multiple Flushes or disinfections may be required. Work out where, fast.

Contact us to discuss Accelerated Recommissioning with Hydrosense with one of our consultants.

Enquire at <https://bit.ly/AcceleratedRecommissioning> or email us at hydrosense@albagaia.com.

#AcceleratedRecommissioning

hydrosense
www.hydrosense-legionella.com